










TAPPER POINT






| Needle Length | Needle Profile | Suture Length | Suture Size | | | | | | | | | | |
|--|------------------|-----------------------------|-------------|-----|-----|----------|----------|-------------|------------|-------------|----------|--|--|
| | | | 8-0 | 7-0 | 6-0 | 5-0 | 4-0 | 3-0 | 2-0 | 0 | 1 | | |
| 16 mm  | ½ Circle | 45 cm | | | | ANW 2303 | | | | | | | |
| | | 70 cm | | | | | ANW 2304 | | | | | | |
| 20 mm  | ½ Circle | 45 cm | | | | | ANW 2304 | | | | | | |
| | | 70 cm | | | | | ANW 2304 | ANW 2437 | | | | | |
| 27 mm  | ¾ Circle Premium | 100 cm | | | | | | | ANW 2511 L | | | | |
| 30 mm  | ½ Circle | 30 cm | | | | | | | ANW 2030 | | | | |
| | | 45 cm | | | | | | | ANW 2317 M | | | | |
| | | 45 cm (5 lengths per foil) | | | | | | | ANW 2508 | | | | |
| | | 45 cm (10 lengths Per foil) | | | | | | | ANW 2503 | | | | |
| | | 70 cm | | | | | | | ANW 2341 | ANW 2342 | | | |
| | | 90 cm | | | | | | | | ANW 2317 | ANW 2338 | | |
| | | | | | | | | | | ANW 2317-3 | | | |
| | | 90 cm Undyed | | | | | | | | | ANW 2337 | | |
| 110 cm | | | | | | | | ANW 2317 VS | | ANW 2827 VS | | | |



| Needle Length | Needle Profile | Suture Length | Suture Size | | | | | | | | |
|--|----------------|----------------------------------|-----------------|-----|-----|-----|---------------|---------------|------------------------------------|---------------|--|
| | | | 6-0 | 5-0 | 4-0 | 3-0 | 2-0 | 0 | 1 | 2 | |
| 40 mm  | 1/2 Circle | 45 cm | | | | | ANW 2345 M | ANW 2346 M | | | |
| | | 45 cm (5 lengths per foil) | | | | | | ANW 2507 | | | |
| | | 45 cm (10lengths Per foil) | | | | | | ANW 2502 | | | |
| | | 70 cm | | | | | | ANW 2546 | | | |
| | | 90 cm | | | | | | ANW 2345 | ANW 2346 ANW 2346-3 | ANW 2360 | |
| | | | 90 cm Undyed | | | | | | ANW 2346 U | | |
| | | 110 cm | | | | | | | ANW 2346 VS ANW 2346 VS-3 | | |
| | | | 120 cm | | | | | | | ANW 2360 L | |

| | | | | | | | | | | | |
|---|-----------------|-----------------|--|--|--|--|--|---------------------|----------------|------------------|---------------|
| <p>40 mm HEAVY</p>  | <p>½ Circle</p> | 20 cm | | | | | | | ANW 2347 SP | | |
| | | 45 cm | | | | | | | | ANW 2347 M | |
| | | 70 cm | | | | | | | | ANW 2359 | |
| | | 90 cm | | | | | | | | ANW 2347 | ANW 2352 |
| | | | | | | | | | | ANW 2347-3 | |
| | | 90 cm Undyed | | | | | | | | ANW 2344 | ANW 2344-3 |
| | | | | | | | | | | ANW 2344-3 | |
| 110 cm | | | | | | | | ANW 2347 VS | ANW 2352VS | | |
| | | | | | | | | ANW 2347 VS-3 | | | |
| 120 cm | | | | | | | | ANW 2350 | | | |
| <p>45 mm HEAVY</p>  | <p>½ Circle</p> | 90 cm | | | | | | | ANW 2347 PB | | |
| | | 120 cm | | | | | | | | ANW 2347 PB L | |
| <p>50 mm</p>  | <p>½ Circle</p> | 75 cm | | | | | | | ANW 9251 | | |
| | | 90 cm | | | | | | | ANW 2346 SP | | |





| | | | | | | | | | | |
|---|-----------------|--------|--|--|--|--|--|--|--|------------|
| <p>50 mm HEAVY</p>  | <p>½ Circle</p> | 75 cm | | | | | | | | ANW 2355 |
| | | 100 cm | | | | | | | | ANW 2448 |
| | | 110 cm | | | | | | | | ANW 2448 L |




TAPER CUT

| Needle Length | Needle Profile | Suture Length | Suture Size | | | | | | | |
|--|----------------|---------------|-------------|-----|-----|-------------|-------------|---|------------|---|
| | | | 6-0 | 5-0 | 4-0 | 3-0 | 2-0 | 0 | 1 | 2 |
| <p>17 mm</p>  | ½ Circle | 45 cm | | | | ANW 2515 | | | | |
| <p>30 mm</p>  | ½ Circle | 45 cm | | | | ANW 2511 HK | | | | |
| | | 70 cm | | | | | ANW 2518 SP | | | |
| <p>30 mm HEAVY</p>  | ½ Circle | 45 cm | | | | | | | ANW 2511 G | |
| <p>36 mm</p>  | ½ Circle | 90 cm (OB-GY) | | | | ANW 2516 | | | ANW 2517 | |
| <p>38 mm</p>  | ½ Circle | 45 cm | | | | | ANW 2511 A | | | |




| | | | | | | | | | | |
|--|----------|---------------|--|--|--|--|--|--|------------|----------|
| 40 mm  | ½ Circle | 45 cm | | | | | | | ANW 1045 | |
| | | 90 cm (OB-GY) | | | | | | | ANW 2518 | ANW 2519 |
| | | 90 cm Undyed | | | | | | | ANW 2518 U | |
| | | 110 cm | | | | | | | ANW 2518 L | |
| 40 mm HEAVY  | ½ Circle | 70 cm | | | | | | | | ANW 2363 |
| | | 90 cm (OB-GY) | | | | | | | | |







CUTTING






| Needle Length | Needle Profile | Suture Length | Suture Size | | | | | | | |
|--|------------------|---------------|-------------|----------|----------|------------|-----|---|---|---|
| | | | 6-0 | 5-0 | 4-0 | 3-0 | 2-0 | 0 | 1 | 2 |
| 16 mm  | ¾ Circle Premium | 45 cm | | ANW 2442 | | ANW 2471 | | | | |
| | | 70 cm Undyed | | | ANW 2443 | | | | | |
| 16 mm  | 5/8 Circle | 45 cm | | | ANW 2465 | | | | | |
| 19 mm  | ¾ Circle Premium | 70 cm Undyed | | | | ANW 2402 S | | | | |
| | | 90 cm Undyed | | | | ANW 2402 | | | | |
| 22 mm  | ¾ Circle | 45 cm | | | | ANW 2401 | | | | |



| | | | | | | | | | | | |
|--|---------------------|-----------------|--|--|--|--|-------------|--|--|--|--|
| 22 mm  | ½ Circle | 90 cm | | | | | ANW 2472 | | | | |
| 25 mm  | ¾ Circle Premium | 90 cm Undyed | | | | | ANW 2328 | | | | |
| 26 mm  | ¾ Circle | 90 cm Undyed | | | | | ANW 2404 | | | | |

REVERSE CUTTING


| Needle Length | Needle Profile | Suture Length | Suture Size | | | | | | | | |
|--|---------------------|-----------------|-------------|-------------|-------------|-----------------------|-------------------|---------------------|---|---|--|
| | | | 6-0 | 5-0 | 4-0 | 3-0 | 2-0 | 0 | 1 | 2 | |
| 12 mm  | ¾ Circle | 45 cm | ANW 2671 | | ANW 2494 | ANW 2495 | | | | | |
| 12 mm  | ¾ Circle Premium | 45 cm | | ANW 2493 | | | | | | | |
| 23 mm  | ½ Circle | 30 cm Undyed | | | | LaNW 2693 SA(S) | | | | | |
| | | 45 cm | | | | | ANW 2510 SB | | | | |
| | | 70 cm | | | | | | ANW 2510 SB/L | | | |
| | | 90 cm Undyed | | | | | ANW 2693 SA | | | | |

| | | | | | | | | | | | |
|--|----------|--------------|--|--|--|----------|------------|---------------|------------|--|--|
| 23 mm HEAVY  | ½ Circle | 35 cm | | | | | | | ANW 2826 | | |
| | | 70 cm | | | | | | | ANW 2826 L | | |
| 26 mm  | ¾ Circle | 70 cm Undyed | | | | ANW 2407 | | | | | |
| 27 mm  | ¾ Circle | 45 cm | | | | | ANW 2694 | ANW 2511 SB | ANW 2827 | | |
| | | 70 cm | | | | | | ANW 2511 SB/L | | | |
| 30 mm  | J Type | 90 cm | | | | | | ANW 2512 | ANW 2828 | | |
| 36 mm  | ½ Circle | 45 cm | | | | | | ANW 2534 M | | | |
| | | | | | | | | ANW 2534 DS | | | |
| | | 90 cm | | | | | | | ANW 2534 | | |
| | | | | | | | | | ANW 2534 T | | |
| 40 mm  | ½ Circle | 45 cm | | | | | ANW 2382 M | | | | |
| | | 90 cm | | | | | ANW 2382 | ANW 2422 | ANW 2478 | | |
| | | 90 cm Undyed | | | | | ANW 2382 U | | | | |

| | | | | | | | | | | |
|--|-----------------|-------------------------|--|--|--|--|-----------------------|----------------------------|------------------------------|---------------------|
| <p>26 mm</p>  | <p>½ Circle</p> | <p>70 cm Undyed</p> | | | | | <p>ANW 2407</p> | | | |
| <p>27 mm</p>  | <p>½ Circle</p> | <p>45 cm</p> | | | | | <p>ANW 2694</p> | <p>ANW 2511 SB</p> | <p>ANW 2827</p> | |
| | | <p>70 cm</p> | | | | | | | <p>ANW 2511 SB/L</p> | |
| <p>30 mm</p>  | <p>J Type</p> | <p>90 cm</p> | | | | | | <p>ANW 2512</p> | <p>ANW 2828</p> | |
| <p>36 mm</p>  | <p>½ Circle</p> | <p>45 cm</p> | | | | | | <p>ANW 2534 M</p> | | |
| | | | | | | | | | <p>ANW 2534 DS</p> | |
| | | <p>90 cm</p> | | | | | | | <p>ANW 2534</p> | |
| | | | | | | | | | <p>ANW 2534 T</p> | |
| <p>40 mm</p>  | <p>½ Circle</p> | <p>45 cm</p> | | | | | <p>ANW 2382 M</p> | | | |
| | | <p>90 cm</p> | | | | | | <p>ANW 2382</p> | <p>ANW 2422</p> | <p>ANW 2478</p> |
| | | <p>90 cm Undyed</p> | | | | | | <p>ANW 2382 U</p> | | |

| Needle Length | Needle Profile | Suture Length | Suture Size | | | | | | | |
|---|----------------|---------------|----------------------|-----|-----|-----|-----|---|---------------|-------------|
| | | | 6-0 | 5-0 | 4-0 | 3-0 | 2-0 | 0 | 1 | 2 |
| 40 mm HEAVY  | 1/2 Circle | 45 cm | | | | | | | ANW 2421 M | |
| | | 90 cm | | | | | | | ANW 2421 | ANW 2351 |
| | | | (1x90cm) (9x45cm) | | | | | | | |
| 50 mm HEAVY  | 1/2 Circle | 90 cm | | | | | | | ANW 2438 | ANW 2353 |
| | | | | | | | | | | |

DOUBLE ARMED

| Needle Length | Needle Profile | Suture Length | Suture Size | | | | | |
|--|---|---------------|-------------|-----|-------------------|-------------------|-------------------|-------------------|
| | | | 4-0 | 3-0 | 2-0 | 0 | 1 | 2 |
| 40 mm Double Armed  | 1/2 Circle Round Bodied, 1/2 Circle TaperCut | 140 cm | | | | ANW 2340 | ANW 2340 | |
| | | 160 cm | | | | | ANW 2160 GF | |
| | | 180 cm | | | ANW 2220 GF | ANW 2080 GF | ANW 2180 GF | ANW 2280 GF |